

Human MSCs Release Multiple EV Populations Containing Mitochondria

CNF Project Number: 2864-20

Principal Investigator(s): Dr. Michelle L. Delco, DVM, Ph.D., DACVS

User(s): Matthew Thomas

Affiliation(s): College of Veterinary Medicine, Department of Clinical Sciences; Cornell University

Primary Source(s) of Research Funding: Harry M. Zweig Fund for Equine Research National Institute of Health-National Institute of Arthritis and Musculoskeletal and Skin Diseases

Contact(s): mld12@cornell.edu, mt826@cornell.edu

Primary CNF Tools Used: Malvern Nano ZS Zetasizer, Malvern NS300 NanoSight

Abstract:

A growing body of evidence supports intracellular mitochondrial (MT) transfer as an important intercellular signaling mechanism. Further, increasing evidence suggests that mesenchymal stromal cells (MSCs) can rescue injured and dysfunctional cells by donating whole mitochondria, and this phenomenon may explain the beneficial effects of therapeutically implanted MSCs. One possible mechanism for MT transfer involves packaging mitochondria into extracellular vesicles (EVs). This would open the possibility of cell-free MT-targeted regenerative therapies. Confirming that this is possible would be an important step toward therapeutic development. As demonstrated here, human MSCs produce EVs containing MT. We have used the CNF to further characterize these 'mitovesicles' and found that there are multiple populations of different sizes, indicating different modes of biogenesis and/or distinct bio-signaling functions.

Summary of Research:

The phenomenon of intercellular mitochondrial transfer, by which mesenchymal stromal cells donate whole mitochondria (MT) to other cell types undergoing MT dysfunction, is a promising avenue for therapeutic intervention in degenerative disease [1]. Mitochondrial donation has been demonstrated in multiple cell types, including neurons and myocytes. It has been shown to improve MT function and prevent apoptosis *in vitro*, as well as improve tissue repair *in vivo* [2-4].

Our lab studies MSC MT donation in the context of orthopedic disease, using *in vitro* chondrocyte cultures and explanted cartilage tissue as models. Using confocal imaging, we have identified several possible modes of MSC-chondrocyte MT transfer, including direct cell-cell contact (e.g., nanotubule-medial filipodial transfer, gap junction-mediated transfer) and what appears to be non-contact transfer, whereby MSCs release mitochondria into the extracellular space, which are then taken up by chondrocytes. We hypothesize that these are MT are packaged inside of extracellular vesicles (EVs) as so-called 'mitovesicles'.

This strategy of loading MT into EVs has precedent in literature. Phinney, et al., showed that MSCs can use mitovesicles, to outsource mitophagy of depolarized MT to

macrophages, boosting bioenergetics for both the donating MSC and the recipient macrophage [5]. Furthermore, Morrison, et al., used cellular staining and flow cytometry to demonstrate that distressed lung epithelial cells can take up MT through EV-mediated transfer, and this ameliorates lung injury *in vivo* [6].

Our goal was to characterize the EVs produced by human MSCs. EVs are an inherently heterogenous population, making specific categorization difficult. However, it is widely recognized that they fall into three size categories: small (15-100 nm, exosomes), medium (150-1000, microvesicles), and large (1 $\mu\text{m}+$, apoptotic bodies).

We isolated EVs from human bone marrow derived MSCs and used dynamic light scattering (DLS) to analyze their size distribution (Figure 1). As expected, we found the three categories supported by previous work [7] (Figure 1).

Next, we used the Malvern NS300 NanoSight to identify which, if any of these categories, contain MT. We stained hMSC EVs with Mitotracker Red, then performed nanoparticle tracking analysis using the Nanosight's 565 nm fluorescent filter. We were able to validate that this method allows us to track exclusively EVs that contain mitochondrial content (Figure 2.)

Further, we found that mitovesicles make up around 20% of the total EV's released and appear to trend slightly larger than the general population (Figure 2.).

The significance of these findings is not yet clear, but likely reflects distinct modes of biogenesis and cargos for different sub-populations of mitovesicles.

Conclusions and Future Steps:

Our work thus far has confirmed our ability to isolate EVs from MSCs and identify mitovesicles within that population. Our next step is to identify and separate the EVs that contain functional and non-functional MT. This will allow us to begin identifying the role that these particles play in intercellular signaling and to further investigate MSC-EV mediated MT transfer.

References:

- [1] Delco, M. L., Bonnevie, E. D., Bonassar, L. J., and Fortier, L. A. Mitochondrial dysfunction is an acute response of articular chondrocytes to mechanical injury. *J. Orthop. Res.* 36, (2018).
- [2] Jiang, D., et al. Mitochondrial transfer of mesenchymal stem cells effectively protects corneal epithelial cells from mitochondrial damage. *Cell Death Dis.* 7, (2016).
- [3] Konari, N., Nagaishi, K., Kikuchi, S., and Fujimiya, M. Mitochondria transfer from mesenchymal stem cells structurally and functionally repairs renal proximal tubular epithelial cells in diabetic nephropathy *in vivo*. *Sci. Rep.* 9, (2019).
- [4] Spees, J. L., Olson, S. D., Whitney, M. J., and Prockop, D. J. Mitochondrial transfer between cells can rescue aerobic respiration. *Proc. Natl. Acad. Sci. U. S. A.* 103, (2006).
- [5] Phinney, D. G., et al. Mesenchymal stem cells use extracellular vesicles to outsource mitophagy and shuttle microRNAs. *Nat. Commun.* 6, (2015).
- [6] Morrison, T. J., et al. Mesenchymal stromal cells modulate macrophages in clinically relevant lung injury models by extracellular vesicle mitochondrial transfer. *Am. J. Respir. Crit. Care Med.* 196, (2017).
- [7] Zhang, X., Hubal, M. J., and Kraus, V. B. Immune cell extracellular vesicles and their mitochondrial content decline with ageing. *Immun. Ageing* 17, (2020).




Figure 1: Dynamic Light Scattering of human MSC-derived extracellular vesicles (EVs) reveals three sub-populations of EVs based on size: small, (SEV; ~ 5-10 nm) medium, (MEV; ~ 100-1000 nm) and large (LEV; 5,000-10,000 nm). N = 3.




Figure 2: NTA supports trend of larger EVs containing MT. EVs were stained with Mitotracker Red and run with either a non-fluorescent (brightfield) filter, or with fluorescence exciting at 565 nm. These groups were compared to non-stained EVs and a double filtered PBS control. Unstained particles were undetectable using the 565 filter. Mitotracker positive EVs seemed to show a trend towards larger size, with nearly all of the smallest population disappearing altogether. N = 1.